



Automated External Defibrillator Weekly Checklist

Months:			
	Date	Name	Cad ID
Week 1			
Week 2			
Week 3			
Week 4			
Supervisor			
Week 5			
Week 6			
Week 7			
Week 8			
Supervisor			
Week 9			
Week 10			
Week 11			
Week 12			
Supervisor			
Single Rescuer Adult, Child & Infant 30 Compressions to 2 Breaths		2 Rescuer Child & Infant 15 Compressions to 2 Breaths (Adults are 30 Compressions)	
Child Pad-Pac 1-8 years old Less than 55 lbs	Adult Pad-Pac At least 8 years old At least 55 lbs	At least 100 Compressions Per Minute	

Weekly checks should include:

- 1) Ensure Ready light is green
- 2) Ensure adult battery/pads are plugged in and within expiration date
- 3) Ensure child battery/pads are available and stored in the case with the AED

Sheriff Patrol AED Serial #: _____



Automated External Defibrillator Weekly Checklist

Months:			
	Date	Name	Cad ID
Week 1			
Week 2			
Week 3			
Week 4			
Supervisor			
Week 5			
Week 6			
Week 7			
Week 8			
Supervisor			
Week 9			
Week 10			
Week 11			
Week 12			
Supervisor			
Single Rescuer Adult, Child & Infant 30 Compressions to 2 Breaths		2 Rescuer Child & Infant 15 Compressions to 2 Breaths (Adults are 30 Compressions)	
Child Pad-Pac 1-8 years old Less than 55 lbs	Adult Pad-Pac At least 8 years old At least 55 lbs	At least 100 Compressions Per Minute	

Weekly checks should include:

- 1) Ensure Ready light is green
- 2) Ensure adult battery/pads are plugged in and within expiration date
- 3) Ensure child battery/pads are available and stored in the case with the AED

Sheriff Patrol AED Serial #: _____